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Table of Contents

Instructions for MA35-UK/MA35-UK-G adapter	
- Handheld-Version	P. 28
- Transmission mode measurement	P. 30
- Handheld-version-set	P. 32
Instructions for using the MA38 and MA38-T	
- Handheld-version	P. 34
- Set for measuring inhomogeneous sample	P. 36
Instructions for using the MA80 adapter	P. 38
Instructions for measuring opaque liquids	P. 40
Instructions for measuring transparent liquids	P. 42
Instructions for measuring translucent liquids	P. 44
Instructions for measuring powder	P. 46

Instructions for measuring colours using the MA35-UK / MA35-UK-G d/8° adapter

 Note: The adapter changes the standard 45°/0° geometry of the sph870/sph900 spectrophotometer into a d/8° geometry according to DIN 5033. This means the samples are illuminated with a diffuse light source created in an integrating sphere with a Barium Sulphate surface.

Attaching the adapter to the probe head:

- Insert the sph870/sph900 probe head into the MA35-UK adapter against the stop and fix by lightly tightening the finger screw.
- Connect the supplied cable to the spectrophotometer and to the adapter.

Important: screwed down the connector only lightly!

The measurement mode of the ColorLite sph870/sph900 automatically switches to probe head MA35-UK mode when the black cable is connected to the device. This will also affect the routine for calibrating.

Calibration

- The black reference and the BAM standard must be unscrewed, as seen below (after calibration to protect both parts screw back together)

Fig.:1 BAM-standard with black reference cap screwed on

Fig.:2 BAM-standard and black reference

- Select in the menu "Calibrate" "BAM-standard" and follow the instructions on the display.
- First calibrate with the black reference, as can be seen on the display. Place the black reference on the adapter (see Fig. 2 – MA35-UK fitted to the and then press the "Enter" key on the device.
- Then place the MA35-UK adapter on the white calibration surface of the BAM-standard and press the "Enter" key again.

 Note: The ceramic surface must be free of any contamination. When cleaning it is recommended to use isopropyl alcohol cleaner. Calibrated before measuring a reference standard, or at least once a day.

Fig.:3 - MA35-UK on black reference

Fig.:4 - MA35-UK and BAM standard

Colour measurement of a sample or standard (reference)

- Samples are measured with the adapter in the same way as with the 45°/0° probe head. The measurements are triggered by pressing the button on the side of the adapter or the "Enter" button on the device.
- For non-homogeneous samples, please use the multiple scan mode with at least three scans. The sample extraction must be repeated between each scan (see the user manual for details).
- The results are automatically displayed after the measurement as an average.
- Use the arrow keys UP/DOWN, to display the results in the different colour systems that are selected in the settings menu.

Instructions for measuring translucent samples in transmission mode

With the help of this accessory translucent liquids and plastic plates are measured in transmission mode. This means the samples are illuminated with a diffused light, generated in an integrating sphere, with a barium sulfate coating. It is important that the samples are always prepared in the same way, so that measurements are reproducible.

Note: Even temperature changes of the sample can affect the reproducibility of results!
Make sure that the arrow points to mark the cells in the direction of the measuring head.

Fitting the standard probe head:

The standard probe head must be inserted into the opening of the accessory using the thumb screw to fix which should be tightened only slightly.

Fitting the probe head adapter MA35-UK:

Fit the probe adapter into the opening provided in the holder and tighten lightly with the thumbscrew. The position should be chosen according to the sample thickness.

Make the connection from the MA35-UK adapter to the sph900/sph870 using the supplied.

Device settings:

The measurement geometry of ColorLite sph870/ sph900 must be set to "CA10-UK". To change the mode, select "Settings" from the main menu. In the window select "Geometry" and press "Enter". Select the next menu "Probe Head". Again here select "CA10-UK" and press "Enter".

After completing translucent liquid measurements you should return the probe head setting to "Auto", so that a connected probe head adapter will automatically be detected (see user manual for details).

Calibration:

For measuring translucent liquids the system is calibrated with distilled water. Fill a cuvette with distilled water and place it in the holder provided. For measuring sheet type materials the system is calibrated with an open beam path.

Note: Make sure that there are no air bubbles in the water.

Select the menu position „Calibrate“ and „Open path / water“ and follow the instructions on the display.

Measuring liquids or solid sheet materials

- After calibration, fill the cuvette with the sample to be measured or alternatively place the material to be measured between probe head and the integrating sphere light source. For inhomogeneous samples use the multiple measurement routine with at least three scans (see also manual). The sample must be moved or newly prepared between each measurement.
- The results are automatically displayed after the measurement as an average.
- Use the arrow keys UP/DOWN, to display the results in the different colour systems that are selected in the settings menu.

Instructions for measuring samples with MA35-UK/ MA35-UK-G d/8° probe head adapter and device holder

 Note: The adapter changes the standard 45°/0° geometry of the sph870/sph900 spectrophotometer into a d/8° geometry according to DIN 5033. This means the samples are illuminated with a diffuse light source created in an integrating sphere with a Barium Sulphate surface.

Attaching the adapter to the probe head:

- Insert the sph870/sph900 probe head into the MA35-UK adapter against the stop and fix by lightly tightening the finger screw.
- Connect the supplied cable to the spectrophotometer and to the adapter.

Important: screwed down the connector only lightly!

The measurement mode of the ColorLite sph870/sph900 automatically switches to probe head MA35-UK mode when the black cable is connected to the device. This will also affect the routine for calibrating.

The MA35-UK measuring head of the handheld instrument is inserted from below into the opening provided in the device holder and hand tighten the locking screw. Then insert the supplied measuring adapter from above on the measuring head. Pay attention to the fact that the adapter is flush with the measuring head.

Calibration

- The black reference and the BAM standard must be unscrewed (seen Fig. 1 and 2). After calibration to protect both parts screw back together.

- Then remove the patch on the measuring head.
- Select in the menu "Calibrate" "BAM-standard" and follow the instructions on the display.
- First calibrate with the black reference, as can be seen on the display. Place the black reference on top of the measuring head (see Fig.3) and then press the "Enter" key on the device.
- Then place the BAM-standard with the white calibration surface on top of the measuring head (see Fig.4) and press the "Enter" key again

Fig. 1 – BAM-standard with black reference cap screwed on

Fig. 3 – black reference on top of the probe head

Fig. 2 – BAM-standard and black reference

Fig. 4 – BAM-standard on top of the probe head

Note: The ceramic surface must be free of any contamination. When cleaning it is recommended to use isopropyl alcohol cleaner. Calibrated before measuring a reference standard, or at least once a day.

Colour measurement of a sample or standard (reference)

- Place the cover on the measuring head. Pay attention to the fact that the measuring aperture is flush with the measuring head.
- Place the sample or standard to be measured on the measuring opening of the aperture
- The measurements are triggered by pressing the button on the side of the adapter or the "Enter" button on the device.
- For non-homogeneous samples, please use the multiple scan mode with at least three scans. The sample extraction must be repeated between each scan (see the user manual for details).
- The results are automatically displayed after the measurement as an average.
- Use the arrow keys UP/DOWN, to display the results in the different colour systems that are selected in the settings menu.

Instructions for measuring inhomogeneous samples with the probe head adapter MA38 and MA38-T

For measuring inhomogeneous samples (for example granules) a Probe head adapter MA38 is needed. This accessory extends the diameter of the measuring spot from 3.5 mm to 38 mm.

It is important to prepare the standard and samples in exactly the same way for each measurement. We recommend defining a method, which can vary depending on the nature of the product. Variations of moisture, density and the optical depth i.e. sample height will affect reproducibility.

Attaching the adapter to the probe head

Insert the standard probe head into the opening of the adapter against the stopper and secure by lightly tightening using the finger screw. Connect the cable of the MA38 to the accessory socket of sph870/sph900.

Device settings

The measurement mode of the ColorLite sph900 automatically switches to probe head MA38 mode when the black cable is connected to the device. This will also affect the routine for calibrating.

Calibration

After changing the mode or when starting to measure the device must be calibrated. How often this needs to be repeated, depends on how accurate the measurements have to be. Temperature changes will cause a slight drift. A calibration will compensate this drift.

- Place the white standard (see Fig.:1) with the white side facing towards the adapter.

Fig.: 1 - MA38 BAM-Standard

- Select from the main menu CALIBRATION and follow instructions on the display.

Note: The white standard tile must be 100% clean before using. We recommend using Isopropyl alcohol cleaner.

- Calibrated before each the measurement of a standard, or at least once a day.

Colour measurement of a sample or standard (reference)

- In the menu Settings adjust the number of "SCANS" to three or more depending on homogeneity of the sample. The probe head have to be put on new on the sample after every measuring.
- Choose between measuring a standard or sample and start with the "Enter" key. The resulting value is displayed automatically as an average.
- Use the arrow keys UP/DOWN, to display the results in the different colour systems that are selected in the settings menu.

- Please note that: by measuring the samples through a glass holder the absolute values will have a systematic error.
- The results can be stored on the device and then later downloaded to a PC using ColorDaTra Basic (optional).
- The Software ColorDaTra Professional (optional) with Online-Mode allows the measurement to be triggered from the PC. The results are then automatically displayed on the PC software.

Instructions for measuring inhomogeneous samples with the probe head adapter MA38 and device holder

For measuring inhomogeneous samples (for example granules) a Probe head adapter MA38 is needed. This accessory extends the diameter of the measuring spot from 3.5 mm to 38 mm.

 **Note:** It is important to prepare the standard and samples in exactly the same way for each measurement. We recommend defining a method, which can vary depending on the nature of the product. Variations of moisture, density and the optical depth i.e. sample height will affect reproducibility.

Attaching the adapter to the probe head

Insert the standard probe head into the opening of the adapter against the stopper and secure by lightly tightening using the finger screw. Connect the cable of the MA38 to the accessory socket of sph870/sph900.

Device settings

The measurement mode of the ColorLite sph870/sph900 automatically switches to probe head MA38 mode when the black cable is connected to the device. This will also affect the routine for calibrating.

Calibration

After changing the mode or when starting to measure the device must be calibrated. How often this needs to be repeated, depends on how accurate the measurements have to be. Temperature changes will cause a slight drift. A calibration will compensate this drift.

- Select from the main menu CALIBRATION and follow instructions on the display.
Note: The white standard tile must be 100% clean before using. We recommend using Isopropyl alcohol cleaner.

- Place the white standard (see Fig.: 1), with the white side facing towards the adapter.

Fig.: 1 - MA38 BAM-Standard

- Calibrated before each the measurement of a standard, or at least once a day.
- Colour measurement of a sample or standard (reference)
- Fill the glass cuvette supplied with sample to be measured. Place the filled cuvette on the middle of the adapter MA38 and cover with the supplied light proof cap.

Fig.: 2 - Light proof cap

- In the menu Settings adjust the number of "SCANS" to three or more depending on homogeneity of the sample. Choose between measuring a standard or sample and start with the "Enter" key. The cuvette must be refilled between each measurement. The resulting value is displayed automatically as an average.
- Use the arrow keys UP/DOWN, to display the results in the different colour systems that are selected in the settings menu.

- Please note that: by measuring the samples through a glass holder the absolute values will have a systematic error.
- The results can be stored on the device and then later downloaded to a PC using ColorDaTra Basic (optional).
- The Software ColorDaTra Professional with Online-Mode (optional) allows the measurement to be triggered from the PC. The results are then automatically displayed on the PC software.

Note: To protect the optical window from soiling place the light proof cap on the adapter when not in use.

Instructions for measuring inhomogeneous samples with the probe head adapter MA80 and device holder

For measuring inhomogeneous samples (for example granules) a Probe head adapter MA80 is needed. This accessory extends the diameter of the measuring spot from 3.5 mm to 80 mm.

 Note: It is important to prepare the standard and samples in exactly the same way for each measurement. We recommend defining a method, which can vary depending on the nature of the product. Variations of moisture, density and the optical depth i.e. sample height will affect reproducibility.

Attaching the adapter to the probe head

Insert the standard probe head into the opening of the adapter against the stopper. The measuring head locks automatically in the measuring head adapter. Connect the cable of the MA80 to the accessory socket of sph870/sph900.

Device settings

The measurement mode of the ColorLite sph870/sph900 automatically switches to the adapter mode when the black cable is connected to the device. This will also affect the routine for calibrating.

 Note: When connecting the MA80 probehead adapter to the device, the MA38 adapter appears on the display as connected accessory. This is not an error because the MA80 has the same system settings as the MA38 measuring head adapter.

Calibration

After changing the mode or when starting to measure the device must be calibrated. How often this needs to be repeated, depends on how accurate the measurements have to be. Temperature changes will cause a slight drift. A calibration will compensate this drift.

For the calibration of the device the large ceramic white standard (see Figure 1) the cover as a black reference (Figure 2) and the cuvette is needed.

Fig.: 1 MA80 BAM-Standard

Fig.: 2 Light proof cap

- Select from the main menu CALIBRATION and follow instructions on the display. The calibration mode MA38 BAM standard is also applicable for calibration with the standard MA80 BAM standard
- Place the light proof cap (black reference) on top of the cuvette and push the enter button.

 Note: It is essential that the cuvette is always placed on the adapter for black calibration.

- Remove the cuvette and place the white standard (see Fig.: 1), with the white side facing towards the adapter and push the enter button.

 Note: The white standard tile must be 100% clean before using. We recommend using Isopropyl alcohol cleaner.

Calibrated before each the measurement of a standard, or at least once a day.

Colour measurement of a sample or standard (reference)

- Fill the glass cuvette supplied with sample to be measured. Place the filled cuvette on the middle of the MA80 adapter and cover with the supplied light proof cap.
- In the menu Settings adjust the number of "SCANS" to three or more depending on homogeneity of the sample. Choose between measuring a standard or sample and start with the "Enter" key. The cuvette must be refilled between each measurement. The resulting value is displayed automatically as an average.
- Use the arrow keys UP/DOWN, to display the results in the different colour systems that are selected in the settings menu.
- Please note that: by measuring the samples through a glass holder the absolute values will have a systematic error.
- The results can be stored on the device and then later downloaded to a PC using ColorDaTra Basic (optional).
- The Software ColorDaTra Professional with Online-Mode (optional) allows the measurement to be triggered from the PC. The results are then automatically displayed on the PC software.

 Note: To protect the optical window from soiling place the light proof cap on the adapter when not in use.

Instructions for measuring opaque liquid samples using the CA10 adapter

Important!

It is important to ensure that the samples are always prepared in the same way, so that measurements are reproducible. Sample temperature can affect the reproducibility! Always fit the cuvettes with the clear side and arrow facing the probe head.

Attaching the sph870/sph900 probe head to the adapter
For the measurement of fluids the probe head must be pushed into the CA10-adapter, so that the sprung inside part is pushed out by 2 – 3 mm and secured with the finger screw. Make sure that no cell is inserted when doing this.

👉 Note: When measuring opaque liquids with the CA10, the device must first be calibrated using the white BAM standard (see calibration below for details).

Device settings:

The measurement geometry of ColorLite sph870/sph900 must be set to 45°/0° (Probe head mode 45°/0° or „Auto“ setting in the measurement geometry – see the user manual). Make sure that the sph870/sph900 is set to calibrate with a working standard in „Settings“ - „System“ the option „use working standard“ must be on (see user manual for details).

Calibration

- To calibrate with the probe head connected to the CA10-holder, a special working standard is included (see Fig.2). This is used to calibrate the device. This working standard must be referenced against the BAM standard before it can be used.

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Abbildung 2 -
BAM-Standard

- Select in the menu „Calibrate“ „BAM-standard“ and follow the instructions on the display. After scanning the BAM the device will call for the working standard to be measured. The probe head must be attached to the CA10 before measuring the working standard.
- The calibration with the working standard should be carried out before the start of a series of measurements, or after a longer break.

Measuring the samples or standards

- After calibration, fill the cell with the sample to be measured. Insert into the adapter into the holder and measure the liquid as a standard or sample. For slightly inhomogeneous samples, please use the multiple scan routine with at least three measurements (see the manual). The sample must be prepared new between each measurement.
- The result of the multiple scans is automatically displayed as an average.

- After the measurement by using the UP/DOWN keys on the spectrophotometer, different colour values are displayed. For example, as absolute values or colour differences between samples of a batch and the reference standard.

Instructions for measuring transparent liquid samples using the CA10-LS adapter

It is important to ensure that the samples are always prepared in the same way, so that measurements are reproducible. Sample temperature can affect the reproducibility! Always fit the cuvettes with the clear side and arrow facing the probe head.

Attaching the sph870/sph900 probe head to the adapter

For the measurement of fluids the probe head must be pushed into the CA10-LS adapter, so that the sprung inside part is pushed out by 2 – 3 mm (see Fig.:1) and secured with the finger screw. Make sure that no cell is inserted when doing this.

Fig.: 1 - CA10 LS with sample

Device settings:

The measurement mode of the ColorLite sph870/ sph900 automatically switches to CA10-LS mode, when the black cable is connected to the device. This will also affect the routine for calibrating.

Calibration:

After changing the mode or when starting to measure the device must be calibrated. How often this needs to be repeated, depends on how accurate the measurements have to be. Temperature changes will cause a slight drift. A calibration will compensate this drift.

- For the Zero point calibration fill a disposable cuvette with distilled water and place the cell in the holder by pulling back the inner part of the probe head which must be connected to the adapter. Important! Make sure the water is air bubble free.
- Select Calibration in the sph870/sph900 menu and then „Water“. Follow the instructions on the display.
- After calibration, fill the cell with the sample to be measured. Insert into the adapter into the holder and measure the liquid as a standard or sample. For slightly inhomogeneous samples, please use the multiple scan routine with at least three measurements (see the manual). The sample must be prepared new between each measurement.
- The result of the multiple scans is automatically displayed as an average.
- After the measurement by using the UP/DOWN keys on the spectrophotometer, different colour values are displayed. For example, as absolute values or colour differences between samples of a batch and the reference standard.

Measuring the samples or standards

- After calibration, fill the cell with the sample to be measured. Insert the adapter into the holder and measure the liquid as a standard or sample. For slightly inhomogeneous samples, please use the multiple scan routine with at least three measurements (see the manual). The sample must be prepared new between each measurement.
- The result of the multiple scans is automatically displayed as an average.
- After the measurement by using the UP/DOWN keys on the spectrophotometer, different colour values are displayed. For example, as absolute values or colour differences between samples of a batch and the reference standard.

Instructions for measuring translucent liquids with the sph870/sph900 and device holder

 Note: It is important to prepare the standard and samples in exactly the same way for each measurement. Differences in the sample temperature can affect the reproducibility.

Attaching the probe head to the stand

Insert the standard probe head of the ColorLite sph870/sph900 with the opening facing upwards into the probe head holder and secure by lightly tightening using the finger screw. The tip of the probe head should be slightly lower than the surrounding ring.

For the measurement of translucent / transparent liquids the ceramic-white standard (spacer) is used as a reflector. The sample is measured through the floor of a distortion-free optical Petri dish.

Device settings:

The measurement geometry of ColorLite sph870/sph900 must be set for this work to "liquids". To change the mode, select "Settings" from the main menu. In the window select "Geometry" and press "Enter". Select the next menu "Probe Head". Again here select "liquids" and press "Enter".

After completing translucent liquid measurements you should return the probe head setting to "Auto", so that a connected probe head adapter will automatically be detected (see user manual for details).

Setting the mode to liquids

Calibration

- To calibrate fill the Petri dish with distilled water. Slowly insert the spacer into the water and place the Petri dish on the middle of the probe head holder. Cover the body with the cap.
- Select in the menu "Calibration" – "Water" and follow the instructions on the device display. Important! Place the reflector precisely in the centre of the Petri dish, make sure no air bubbles are trapped between the white reflector and the dish floor (see Fig.: 2).


Fig.:2 Position of the reflector in the Petri dish

Measuring a sample or standard

- After calibration, fill the Petri dish with the sample and insert the white reflector, place on the middle of the probe holder. Cover the dish with inserted reflector and sample with the light proof cap. Measure the liquid as a standard or sample.

- For measuring inhomogeneous samples adjust the number of "Scans" to three or more in the "Settings" menu. The Petri dish must be refilled between each measurement. The resulting value is displayed automatically as an average (see the user manual for details).
- Use the arrow keys UP/DOWN, to display the results in the different colour systems that are selected in the settings menu.

 Note: By measuring the samples through a glass holder the absolute values will have a systematic error.

Instructions for measuring powder with the sph870/sph900

👉 Note: The sample must be prepared each in exactly the same way to achieve a good reproducibility. Moisture and the density of the different samples will affect the reproducibility!

Device settings

Powders are normally measured through the glass floor of the optical cuvette as a solid. The mode of the ColorLite sph870/sph900 should be switched to $45^\circ/0^\circ$ or $d/8^\circ$ accordingly.

👉 Note: The glass floor will deviate the light, which results in a systematic error. A compensation routine is available by ColorLite if needed.

Calibration

- For the daily calibration the large ceramic white standard is used. This must be previously measured against the BAM standard. Select in the menu „Calibrate“ „BAM-standard“ and follow the instructions on the display.
- When asked to measure the working standard first connect the probe head to the holder, fix so that it is not protruding and is slightly lower than holder. Place the large ceramic working standard on the probe head (Fig. 1-sample holder head and large ceramic-white-standard) and press the „ENTER“ button on the device.
- After this the device is calibrated directly using the large ceramic working standard – Select in the menu calibration „Working Standard“.


Measurement

- The sample is measured through the distortion-free floor of the optical cuvette
- Fill the powder in the cuvette approximately up to a height of 10 mm. To spread the powder evenly, lightly tap the filled cuvette on a hard surface. Then place on the middle of the holder and cover with the light proof cover.
- Slightly inhomogeneous samples should be measured using the multiple scanning mode. With the sample being prepared new between each scan. (see the user manual for more details).

Note: A special set (MA38 see P.63) is available from ColorLite for measuring inhomogeneous samples.

If using the multiple scan mode is used the results are displayed as an average value. Use the arrow keys UP/DOWN, to display the results in the different colour systems that are selected in the settings menu. For example the absolute values or colour differences between the standard and sample Also it is possible to visualise the results with the PASS/FAIL output.

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ColorLite GmbH